



CELLMED METHOD STEP 1 – OPENING UP THE SKIN

## SKIN BOOSTER SENSITIVE pH 3.0



Sensitive skin reacts quickly by becoming irritated and turning red. It's often the case here that sun and cosmetics are not well tolerated, and stress can also have a rapid negative effect on skin complexion and condition. SKIN BOOSTER SENSITIVE pH 3.0 Acid Treatment was developed especially for sensitive skin. Aligned with the needs of sensitive skin, it gently loosens dead skin cells and supports skin regeneration. It refines skin structure and increases the skin's moisture-retention capability.

PROFESSIONAL

### AREA OF APPLICATION

- Sensitive skin prone to irritation and reddening

### EFFECT

- Removes dead skin cells and stimulates skin cell renewal
- Has a moisture-retention effect
- Promotes skin cell renewal
- Prepares the skin for subsequent skincare products
- When used regularly, it boosts the resilience of sensitive skin
- Helps reduce the intensity of expression lines

### MAIN INGREDIENTS

AHA acid (lactic acid)

This acid has a larger molecular structure than glycolic acid, which means it is unable to penetrate that deeply into the skin and is thus perfect for use on sensitive skin. It removes dead skin cells, stimulates skin cell renewal, and optimizes the skin's moisture-retention capability.

PHA acid (lactobionic acid)

gently removes dead skin cells, optimizes the skin's moisture-retention capability, and counteracts oxidative processes

Hexapeptides

have a muscle-relaxing effect and reduce the intensity of expression lines

Bisabolol

Mild and soothing properties

### APPLICATION

Make sure to follow the treatment plan for the BIODROGA MD Acid Treatment and pay attention to the contraindications.

### PRODUCT LINKS

- SKIN BOOSTER Acid Neutralizer
- HYPER SENSITIVE line

Product rating  
[www.kosmetikanalyse.org](http://www.kosmetikanalyse.org)  
"Good"